Phonak hearing aids

Our hearing aids bring you the best possible solutions for your hearing needs. **Hearing Aids Latest**



Phonak Hearing Aids

Virto B-Titanium

By incorporating a hightech material for Virto B-Titanium, we've created hearing aids that reflect our passion for innovation and ingenious engineering design. This new Phonak IIC combines performance, functionality, and aesthetics. With a significantly thinner shell and 60% smaller electronic components, we are able to reduce the size of Virto B-Titanium by up to 26% compared to its predecessor. With more space available inside the hearing aid, a larger vent can be integrated, reducing occlusion for more listening comfort. Alternatively more space allows for the integration of a superpower receiver, resulting in an impressive 60% fit rate for IIC devices.





Audéo B

Phonak Audéo B hearing aids were created to simplify your life and make hearing effortless. What makes Audéo B so special is its unique operating system called AutoSense OS. It senses your surroundings and adjusts your hearing aids every step of the way. All you have to do is turn them on, and the rest is automatic. And with the rechargeable Audéo B-R, you can now enjoy 24 hours* of hearing with one simple charge.

Andéo B-A

It doesn't matter where you are or what you're doing, Phonak rechargeable hearing aids were created to simplify your life. Finally, a chance to free yourself from the hassles of disposable batteries with a groundbreaking built-in lithium-ion battery technology. And because it's a phonak, you can always be assured of superior hearing performance of superior hearing performance

A-d obbuA





Phonak NaídaTM V

Phonak Naída V utilizes the Venture platform with a host of new groundbreaking dedicated features to improve audibility for people with severe to profound hearing loss. With more power than ever before, more people can benefit from a smaller device



Bolero V

The Finest in Performance and Reliability

Phonak Bolero V Behind-the-Ear (BTE) hearing aids feature our latest high-performance technology and are designed to provide you a seamless listening experience, even in the most challenging environments. In an appealing and robust design, Bolero V is available in four performance levels and offers a wide choice of colors to either blend in or stand out from your hair or skin tone.



Bolero V

The Finest in Performance and Reliability

Phonak Bolero V Behind-the-Ear (BTE) hearing aids feature our latest high-performance technology and are designed to provide you a seamless listening experience, even in the most challenging environments. In an appealing and robust design, Bolero V is available in four performance levels and offers a wide choice of colors to either blend in or stand out from your hair or skin tone.

Wireless Accessories



Roger Pen is a versatile wireless microphone that reduces loud background noise and transfers the clear sound of the voice you want to hear directly to your hearing aids.



Roger EasyPen is designed for you to interact and understand in challenging listening situations. Packed in a stylish design, it features a fully automatic wireless microphone which is easy to use in one-on-one or group conversations.



Roger Clip-On Mic is a discreet shirt-worn wireless microphone that is worn by your partner for one-to-one communication in loud noise and over distance.



DECT - looks like a regular cordless phone, but is actually much more powerful. It transmits sound to both your hearing aids simultaneously, reducing noise and maximizing understanding.



EasyCall - Phonak EasyCall connects any Phonak wireless hearing aid wirelessly to your cell phone. It is permanently attached and connected to the phone and streams the conversation directly to both hearing aids.



The ComPilot Air is a small and stylish streamer that offers excellent stereo sound quality while connecting Phonak hearing aids via Bluetooth to a variety of audio sources.



RemoteMic is used together with a Phonak streamer to understand a speaker over distance. It is attached to your conversation partner's clothing and transmits their voice to your hearing aids.



TVLink - turns your Phonak hearing aids into a wireless TV headset with stereo sound quality.

